| | | |
|---|--|---------------------|
|  | 记录编号 FN: SMP-RD-A0002-R06 | 版本号 Ver: 2.1 |
| | 归档 Filing: filed immediately, permanently kept | 生效日期 Ed: 2023.01.12 |

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Medium for Animal-Cells

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Maxgrow 201 is a chemically-defined powder medium designed especially for the growth of Chinese Hamster Ovary (CHO) cells and production of recombinant protein in CHO cells. The medium is supplied without L-glutamine for use in Glutamine Synthetase, or GS System™. For other applications, add low levels of L-glutamine or GlutaMAX™ prior to use. The recommended glucose concentration to add into Maxgrow 201 Medium is 8g/L under typical culture conditions and won't cause potential lactic acid accumulation problem.

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This product is for research or for further manufacturing use. This product is not intended for human or therapeutic use.

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Suitable for Maxgrow 201 powder packaging (take 1L as an example)

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1. Prepare ultrapure water with 80% of the volume of the solution (room temperature 20~30℃);

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Recommended cell passage density: $0.3\sim 0.5\times 10^6$ cells/mL, when the cell density reaches $1.5\sim 3.5\times 10^6$ cells/mL, the passage operation is carried out.

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In many cases, CHO cells such as CHO-K1, CHO-K1 SV, CHO-S may be subcultured from a serum-free medium or chemically defined medium directly into Maxgrow 201.

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1. Dispense Maxgrow 201 medium into an appropriate culture vessel and equilibrate to 37° C and 5% CO₂.

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2. Passage CHO cells from serum-free culture or other chemically defined media into Maxgrow 201 at $3 - 5 \times 10^5$ viable cells/mL. It is important that cells should be in the logarithmic phase of growth with at least 90% viability before passaging.

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3. Incubate cultures at 37°C and $5\% \text{CO}_2$ until the viable cell density greater than 1×10^6 cells/mL, preferably $1.5 - 3.5 \times 10^6$ cells/mL.

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4. Subculture into Maxgrow 201 medium at $3 - 5 \times 10^5$ cell/mL seeding density.

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5. Maintain cells in Maxgrow 201 for at least three passages to allow complete adaptation and assure optimum performance.

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6. For the culture required L-Glutamine supplementation, titrate L-Glutamine concentration from 1mM to 4mM for the best results.

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Freezing solution: 45% Maxgrow 201 + 45% condition medium + 10% DMSO. No proteins or other additions are necessary.

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Collect the cells in a centrifuge tube, centrifuge at 800-1000rpm for 5 minutes and discard the supernatant;

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Resuspend the cells in the freezing solution, store $1.0 \times 10^7 - 1.5 \times 10^7$ cells in each cryovial, put them in a freezing box, and transfer them to a liquid nitrogen tank for permanent storage after overnight at -80°C .

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① Prepare 37°C warm water, take the cells from liquid nitrogen and put them in warm water immediately to quickly thaw the cells;

| | | |
|---|--|---------------------|
|  | 记录编号 FN: SMP-RD-A0002-R06 | 版本号 Ver: 2.1 |
| | 归档 Filing: filed immediately, permanently kept | 生效日期 Ed: 2023.01.12 |

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② Directly resuspend the cells in one cryopreservation tube in pre-warmed Maxgrow 201, control the cell density between $3E5$ - $5E5$, and transfer to a 125 ml shake flask for culture;

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④ After 2-3 days of culture, count and passage the cells.

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Cells were passaged at the density of $3E5$ - $5E5$, counted every 2-3 days, and passaged. For the first 3 passages, the volume was unchanged to restore cell viability. The standard of normal seed status: the vigor is more than 90%, the cell shape is regular and round, and the growth doubling time is normal.

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Maxgrow 201 powder packaging: store in dark at $2\sim 8^{\circ}\text{C}$, and the retest period is 24 months.

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